



# AVID FriendChips (SUDS)

NOW Available in: 15ga/8mm x 1.4mm or 12ga/12mm x 2.1mm

## FriendChip Mini -15ga (25 ct)

AVID 20881 (Avid Code), AVID23881 (ISO)

## FriendChip - 12ga (25 ct)

AVID 2028 (Avid Code), AVID2128 (FECAVA), AVID2328 (ISO)

## About the AVID FriendChip

- Operating Temperature: -25 to 80C
- Frequency: Available in 125 kHz or 134.2 kHz
- ISO Compatible (FDX-B)
- Housing Material: Bioglass
- EMC: Not affected by X-Ray or MRI
- Passive Power: No batteries, reader activated
- Read Distance: Depends on reader 8-20 cms
- Life Expectancy: 1,000,000 (+) read activations

Size: 15ga/8mm x 1.4mm or 12ga/12mm x 2.1mm

Avid's 12ga or 15ga SUDS delivery systems are preloaded and sterilized, and ready for use.



This (SUDS) - Single Use Delivery System incorporates parts that are custom made for AVID by Monoject. Each unit features corresponding identification tag and four peel and stick labels for easy, accurate documentation and record keeping.

Call Now Toll-Free **(800) 336-2843** Sales Option 3  
<http://AvidID.com>