



AVID FriendChips (SUDS)

NOW Available in: 15ga/10mm x 1.4mm or
12ga/12mm x 2.1mm

FriendChip Mini -15ga (25 ct)

AVID 20881 (Avid Code), AVID23881 (ISO)

FriendChip - 12ga (25 ct)

AVID 2028 (Avid Code), AVID2128 (FECAVA), AVID2328 (ISO)

About the AVID FriendChip

Operating Temperature: -25 to 85°C

Frequency: Available in 125 kHz or 134.2 kHz

ISO Compatible (FDX-B)

Housing Material: Bioglass

EMC: Not affected by X-Ray or MRI

Passive Power: No batteries, reader activated

Read Distance: Depends on reader 8-20 cm

Life Expectancy: 1,000,000 (+) read activations

Size: 15ga/10mm x 1.4mm or
12ga/12mm x 2.1mm

Avid's 12ga or 15ga SUDS delivery systems
are preloaded and sterilized, and ready for use.



This (SUDS) - Single Use Delivery System incorporates parts that are custom made for AVID by Monoject. Each unit features corresponding identification tag and four peel and stick labels for easy, accurate documentation and record keeping.

Call Now Toll-Free **(800) 336-2843** Sales Option 3
<http://AvidID.com>