



# AVID FriendChips (SUDS)

NOW Available in: 15ga/8mm x 1.4mm  
15ga/10mm x 1.41mm  
12ga/12mm x 2.1mm

## FriendChip Mini - 15ga (25 ct)

AVID20881 (Avid Code), AVID23881 (ISO)

## FriendChip - 12ga (25 ct)

AVID2028 (Avid Code), AVID2128 (FECAVA), AVID2328 (ISO)

## About the AVID FriendChip

- Operating Temperature: -25 to 85°C
- Frequency: Available in 125 kHz or 134.2 kHz
- ISO Compatible (FDX-B)
- Housing Material: Bioglass
- EMC: Not affected by X-Ray or MRI
- Passive Power: No batteries, reader activated
- Read Distance: Depends on reader 8-20 cm
- Life Expectancy: 1,000,000 (+) read activations
- FriendChip Size: 12ga/12mm x 2.1mm
- FriendChip Mini Sizes: 15ga/8mm x 1.4mm for AVID20881 (Avid Code)  
15ga/10mm x 1.41mm for AVID23881 (ISO)

Avid's 12ga or 15ga SUDS delivery systems are preloaded and sterilized, and ready for use.



This (SUDS) - Single Use Delivery System incorporates parts that are custom made for AVID by Monoject. Each unit features corresponding identification tag and four peel and stick labels for easy, accurate documentation and record keeping.

Call Now Toll-Free **(800) 336-2843** Sales Option 3  
<http://AvidID.com>